

Taps, saddle valve and tubing

Mini-taps dispense purified water at the point-of-use. A simple, stylish and elegant design that looks great in any application and especially a drinking water outlet in the kitchen.

Attaching the tap to the existing water pipe is made simple with the use of a saddle valve to pierce and clamp the pipe and redirect water via the filter to the tap. Tubing is used in this application to connect the filter to the tap and many water systems require tubing to connect equipment. Tubing is available in ¼", 3/8" and ½" sizes in either full rolls or cut to metre lengths.



Mini Tap - MT1



Saddle valve - SV1



Tubing

Product information					
Product Code	Description	Box Quantity	Roll Length (m)	Box Weight (kg)	Roll Weight (kg)
MT1	Mini Tap - Chrome	25	-	8	-
SV1	Saddle tapping valve with tubing	6	-	2	-
¼" Tube	Blue/Red/White/Black/Natural 150m	-	150	-	2
¾" Tube	Blue/Red/White/Black/Natural 150m	-	150	-	2
½" Tube	Blue/Red/White/Black/Natural 75m	-	75	-	2

TDS Meter

TDS meters are used to measure the levels of Total Dissolved Solids (TDS) in a given volume of water, an essential factor affecting filter performances and required in the correct set-up for reverse osmosis systems.



TDS-3 – A simple, hand held, slimline unit, which comes with a carry case and has the additional feature of being a thermometer.

TDS-DM2 – A dual inline monitor allowing the measurement of two different water lines at any time. This highly accurate, digitally calibrated model comes with diagnostic error messaging, shielded sensor cables, water resistant casing and is ideal for industrial applications.

Product information			
Product Code	Description	Box Quantity	Weight (kg)
TDS-3	Hand held digital TDS meter	10	1
TDS-DM2	Dual TDS meter c/w ¼" JG connections	6	1

Countdown Timer

The CT1 from SPECTRUM is a useful accessory to help remind users when a changeout is needed. Suitable for a variety of applications the CT1 can be used to set reminders for cartridge and UV lamp changeouts to ensure systems run efficiently and effectively.



Product information			
Product Code	Description	Box Quantity	Weight (kg)
CT1	Countdown timer	20	1